

## DIRECTIONS FOR REMOVAL:

Within 24 hours, soak GlykoSL<sup>®</sup> with water and gently remove.


Store at controlled room temperature,  
15 °C to 30 °C (59 °F to 86 °F).

Sterilized with radiation.

Do not resterilize or reuse.

Single use only.

Not Made with Natural Rubber Latex



MPS 500-121 rev 03 05/2018

# GlykoSL<sup>®</sup>

REF 400-24


REF 400-26

## MPT MARINE POLYMER TECHNOLOGIES

Innovative Medical Solutions

159 Lorum St.  
Tewksbury, MA 01876-1716  
1-888-666-2560

Made in the U.S.A.



**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

**DESCRIPTION:**

GlykoSL® is a conformable, soft and absorbent calcium alginate dressing.



**INDICATIONS:**

GlykoSL® is indicated for use in the local management of bleeding wounds such as vascular access sites and percutaneous catheters or tubes.

**PRECAUTIONS:**

There are no known contraindications to the use of GlykoSL®.  
Longer compression may be necessary for hypertensive or obese patients.

**DIRECTIONS FOR APPLICATION:**

1. Open GlykoSL® pouch into sterile field.
  2. Hold pressure proximal to the puncture site.
  3. Place GlykoSL® over puncture site.
  4. Apply firm compression over the puncture site and release the proximal pressure after 3 to 4 minutes.
  5. Hold pressure on puncture site until hemostasis is achieved, then slowly release.
  6. Place dry gauze over GlykoSL® and cover with a transparent dressing.
- 