




EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



Syvek is a registered trademark.

US Patent Nos.

7,041,657, 7,115,588, 7,931,637, 8,152,750,  
8,481,512, 8,834,516, 8,835,408, 8,859,528,  
8,992,453

MPS 500-62 REV-06 5/2018

**Syvek**Excel®  
PATCH

REF 400-16

**MPT** MARINE POLYMER  
TECHNOLOGIES

Innovative Medical Solutions



159 Lorum St.  
Tewksbury, MA 01876-1716  
1-888-666-2560

Made in USA

## DESCRIPTION

SyvekExcel® Patch is a soft, white, sterile, lyophilized pad of a cellulosic polymer (poly-N-acetylglucosamine) isolated from microalgae, and attached to medical grade foam backing.

## INDICATIONS

SyvekExcel® Patch is intended for use under the direction of a healthcare professional.

SyvekExcel® Patch promotes the rapid control of bleeding in patients following hemodialysis and in patients on anticoagulation therapy.

SyvekExcel® Patch is indicated for use in the local management of bleeding wounds such as vascular access site, percutaneous catheters or tubes and surgical debridement.

## PRECAUTIONS

There are no known contraindications to the use of SyvekExcel® Patch. Longer compression may be necessary for hypertensive or obese patients.

### **Not Made with Natural Rubber Latex**

Store at controlled room temperature

15°C to 30°C (59°F to 86°F).

Sterilized with radiation.

Do not resterilize and/or reuse.

Single use only.

## SyvekExcel® Patch APPLICATION

1. Open SyvekExcel® Patch tray into sterile field.
2. Hold occlusive pressure proximal to the puncture site and remove the sheath.
3. Allow 1-2 drops of blood from site.
4. Immediately place SyvekExcel® Patch (beige side AWAY from skin) over the puncture site.
5. Apply firm pressure over the puncture site.
6. Slowly decrease amount of proximal pressure after 4 to 6 minutes (3 - 4 minutes for diagnostic procedures) while continuing to apply firm compression over the puncture site.
7. Hold firm pressure on the puncture site until hemostasis is achieved, then slowly release.
8. Assess the site without disturbing the SyvekExcel® Patch, then cover with a transparent adhesive dressing.

## REMOVE SyvekExcel® Patch after 24 Hours:

1. Remove the transparent adhesive dressing.
2. Soak SyvekExcel® Patch with water.
3. Gently peel off SyvekExcel® Patch.

**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.